



Quick Reference Guide

DHAtoD86 1.1.6

www.StillPeaks.com



StillPeaks DHAtoD86 software consists of 3 apps that read peak integrator data from Chromatographic Data Systems (CDS) and calculate a D86 report of your sample:



Virtual GC is a computer model of a Dimethylsilicone coated capillary gas chromatographic column. Kovats index databases and n-Paraffin retention times are calculated from GC configuration input and an isothermal Kovats index database. The GC configuration should match the actual GC configuration. Column aging is accounted for with film thickness. Copy the Virtual GC window to your lab journal.

DHA Detailed Hydrocarbon Analysis software reads the GC integration file. Peaks are identified by shortlist, pattern and Kovats index database. From DHA, Virtual GC and D86 can be started.

D86 The Fugacity-Film model uses the DHA output of the sample to calculate the distillation curve of the sample. The D86 report is shown on screen and written to the sample folder. To keep the same numbers, you can enter D86 distillation data and offsets are calculated. D86 set points and off sets are managed per sample type.


Content

1.DHA Detailed Hydrocarbon Analysis	2
2.Virtual GC: <i>n</i> -Paraffin Times & Kovats Index	3
3.DHA Peak Identification.....	5
4.D86 Sample Type Setpoints.....	6
5.D86 Property Units and Thermometer Shift	8
6.D86 Calculation and Setpoint Adjustents	8
7.D86 Physical Data Input	13


Quick Reference Guide

1.DHA Detailed Hydrocarbon Analysis

The DHA opening screen is empty:




The Types tabsheet allows you to select your sample type. The typical shortlist of big peaks is displayed:



Use *i*-Paraffins for Alkylates and Isomerates.


Sample type selection also changes the Kovats index database for peak identification:



The Kovats index database does not contain the *n*-Paraffins (indices follow from carbon number I=100c).

Quick Reference Guide


Use directory tree at the left side to select your ChemStation or OpenLab analysis sequence. The chromatogram at the right side displays the integrator results of a selected sample in the samples list:



The white slice of the small circle at the left bottom of the chromatogram indicates almost 50% of the sample peaks are not identified.

2.Virtual GC: *n*-Paraffin Times & Kovats Index

Click the nParaffins tab sheet of the left hand side tab lineup:




n-Paraffin retention times are displayed, along with a Virtual GC button and a Read Alkanes button.


The right side tab sheet lineup shows the measured GC peak times and peak Area of sample RFG 0210 located in the path displayed in the window title. The peak identity result is shown in the 4th column.

The measured first peak time in the GC table is before CH₄ time, hence the *n*-paraffin times needs to be adjusted. Please lookup the GC method used to acquire the GC data. The Virtual GC button enables you to adjust the methane time at 35°C to the 0.35 minute setpoint:

Quick Reference Guide




Select the Save Method & Kovats index data menu option of Virtual GC. Times and databases are filed. Select the Print option of the Virtual GC menu to keep track of the n-Paraffin calculations. The n-Paraffin times and databases are applied in DHA only after clicking the Read Alkanes button:



Quick Reference Guide

3.DHA Peak Identification

The Peak ID pie chart has a red part for the big peaks, a yellow part for the pattern recognition, a blue part for the small peaks with Kovats index and a white part for the unidentified peaks. In this sample 13 peaks have Area% > 2, corresponding 35% red part of the pie chart:




The left hand tab sheet Print is selected, where the hard coded I8 pattern is selected. This pattern is recognized in this sample and marked with I8. The yellow part of the ID pie indicates that about 20% of the peaks in this sample are identified by peak recognition of the Iso-octanes (I8), Xylenes and Ethylbenzene (A8, A9).

NOTE The D86 curve button appears only if less than 50% of the peaks is unknown.

Quick Reference Guide


4.D86 Sample Type Setpoints

Clicking the D86 button opens the D86 page, showing the sample and available D86 data:



The sample has been recognized as Summer Spec gasoline and the ASTM D86 cross check data are displayed at the left side. The setpoints and offsets for this summer gasoline are found in the D 86 Offset tabsheet:

ASTA can be switched off by unchecking the Automatic Sample Type Assignment checkbox:



Quick Reference Guide

Chromatogram | D86 | D 86 Offset | Units

Automatic Sample Type Assignment

Stream Selector	from sample	Tsample[°C]	Theater[°C]
<input type="radio"/> Gas		0	0
<input type="radio"/> Alkylate	Alkylate	15	160
<input type="radio"/> Isomerate	Isomerate	15	500
<input type="radio"/> Light Naphtha	LSR	15	350
<input type="radio"/> Medium Naphtha	RFG 0402	10	400
<input type="radio"/> Heavy Naphtha	HSR	15	400
<input type="radio"/> Light FCC	LFCC	5	400
<input type="radio"/> Medium FCC	MFCC	15	500
<input type="radio"/> Heavy FCC	RFG 0402	25	570
<input type="radio"/> Reformate	MFCC	5	500
<input type="radio"/> Gasoline Winter	RFG 0210	15	450
<input checked="" type="radio"/> Gasoline Summer	RFG 0303	5	350
<input type="radio"/> RBOB	RFG 0210	15	500

RVP sampling temperature 25.0 °C

D86 rate [mL/min]	4.5
Delay?[mL]	4.0
Psystem [Pa]	101325
Print D86 report	N

The ASTA check will stay checked or unchecked.

Now with ASTA unchecked, all sample type specific setpoints can be selected allowing you to run a gasoline with reformate setpoints.

The setpoint selection and the value of the setpoints is retained in the D86 software.

Tsample is the Sample Temperature, Theater is the heater temperature before IBP.

The offset shown at the right side can be applied to the final D86 report by checking the Offset checkbox:

c:\hpchem\1\data\b\007F0701.D\DHA_RFG 0210.csv

File Print Logbook

1.0

Still Peaks

Offset

RVP D 323

Sample [kPa]

RFG 0210 49.66 s

Evaporated Cut V%	Td86 °C	StDev °C
IBP	44.6	2.5
5	56.4	0.8
10	59.2	0.6
20	62.8	0.5
30	65.9	0.5
40	74.3	1.9
50	98.1	0.9
60	109.0	1.0
70	121.5	1.0
80	139.1	0.9
90	163.3	1.1
95	179.9	1.8
FBP	206.1	2.6

Type: Gasoline Summer

Tsample (5°C)	Loss %	-
Theater (350°C)	Recovery %	-
Power	Residue %	-
Time before IBP	Time rec 5%	-
Rate (4.5mL/min)	Time last 5%	-

Chromatogram | D86 | D 86 Offset | Units

Automatic Sample Type Assignment

Stream Selector	from sample	Tsample[°C]	Theater[°C]
<input type="radio"/> Gas		0	0
<input type="radio"/> Alkylate	Alkylate	15	160
<input type="radio"/> Isomerate	Isomerate	15	500
<input type="radio"/> Light Naphtha	LSR	15	350
<input type="radio"/> Medium Naphtha	RFG 0402	10	400
<input type="radio"/> Heavy Naphtha	HSR	15	400
<input type="radio"/> Light FCC	LFCC	5	400
<input type="radio"/> Medium FCC	MFCC	15	500
<input type="radio"/> Heavy FCC	RFG 0402	25	570
<input type="radio"/> Reformate	MFCC	5	500
<input type="radio"/> Gasoline Winter	RFG 0210	15	450
<input checked="" type="radio"/> Gasoline Summer	RFG 0303	5	350
<input type="radio"/> RBOB	RFG 0210	15	500

RVP sampling temperature 25.0 °C

D86 rate [mL/min]	4.5
Delay?[mL]	4.0
Psystem [Pa]	101325
Print D86 report	N

With this check, the offsets are applied to the final D86 calculation. Offset is found by entering ASTM D 86 data and can be applied to tune DHAtoD86 even closer to your D 86 numbers.

Quick Reference Guide

5.D86 Property Units and Thermometer Shift

The units of the properties P, T and density can be selected in the Units tabsheet.

Eö is a dimensionless number, related to droplet size. The bigger the droplet at the virtual thermometer tip, the lower the temperature.

The Eö number shifts the calculated D86 curve down. Any changes are saved and applied ever after. Droplet size can be restored to default values by hitting the default values button.


6.D86 Calculation and Setpoint Adjustents

The Eö number from chapter 5 is much too high, shifting the calculated D86 curve down:


%	Evap* [°C]	Rec* [°C]
IBP	37.4	37.4
5	42.6	44.5
10	48.4	49.8
20	58.4	59.0
30	63.2	63.5
40	65.8	66.0
50	69.8	71.6
60	103.5	105.2
70	116.5	118.0
80	134.6	136.9
90	159.7	163.7
95	178.5	185.6
FBP	204.8	204.8

So we hit the default values button (chapter 5, Eö = 1) and re-run the D86:

Quick Reference Guide




Maybe Eö = 1.666 is better:




Now we have 10 temperatures within the $-2 < z < 2$ range. The IBP is overestimated, so we take an initial heater temperature of 400°C:

Quick Reference Guide



The repair of the 10% may be found by adding another 25°C to Theater:



Next we change the Delay setpoint and the Eo number in a last effort to get closer to the ASTM D 86 data:

Quick Reference Guide

1.0

Still Peaks 


Offset
Sample: RFG 0210 [kPa] 49.66 s

Evaporated Cut V%	Td86 °C	StDev °C
IBP	44.6	2.5
5	56.4	0.8
10	59.2	0.6
20	62.8	0.5
30	65.9	0.5
40	74.3	1.9
50	98.1	0.9
60	109.0	1.0
70	121.5	1.0
80	139.1	0.9
90	163.3	1.1
95	179.9	1.8
FBP	206.1	2.6

Type: Gasoline Summer
 Tsample (5°C) 211 Loss % -
 Theater (425°C) 425 Recovery % 97.05
 Power (14.0W) 14.0W Residue % 1.09
 Time before IBP 8.43 min Time rec 5% 54s
 Revap (4.5mL/min) 3.06 Time last 5% 72s


Chromatogram D86 D 86 Offset Units

RFG 0210
DHA_RFG 0210



Evap	GLC[°C]	D86[°C]	G-D[°C]	z (12ok)
IBP	43.6	44.6	-1.0	-0.4
5	56.7	56.4	0.3	0.4
10	60.6	59.2	1.5	2.5
20	63.6	62.8	0.8	1.5
30	65.5	65.9	-0.3	-0.7
40	70.7	74.3	-3.6	-1.9
50	99.9	98.1	1.8	2.0
60	110.2	109.0	1.2	1.3
70	121.5	121.5	-0.1	-0.1
80	137.6	139.1	-1.5	-1.6
90	162.4	163.3	-0.9	-0.8
95	182.0	179.9	2.1	1.1
FBP	211.1	206.1	5.0	1.9

Refresh Offset



15.09.2004 15:47:09 c:\hp\chem\1\data\b\007f0701.D\DHA_RFG 0210.csv


To repair the last differences, offset is refreshed by clicking the Refresh Offset, and we the result is:

Chromatogram D86 D 86 Offset Units

Automatic Sample Type Assignment


Stream Selector	from sample	Tsample[°C]	Theater[°C]
<input type="radio"/> Gas		0	0
<input type="radio"/> Alkylate	Alkylate	15	160
<input type="radio"/> Isomerate	Isomerate	15	500
<input type="radio"/> Light Naphtha	LSR	15	350
<input type="radio"/> Medium Naphtha	RFG 0402	10	400
<input type="radio"/> Heavy Naphtha	HSR	15	400
<input type="radio"/> Light FCC	LFCC	5	400
<input type="radio"/> Medium FCC	MFCC	15	500
<input type="radio"/> Heavy FCC	RFG 0402	25	570
<input type="radio"/> Reformate	MFCC	5	500
<input type="radio"/> Gasoline Winter	RFG 0210	15	450
<input checked="" type="radio"/> Gasoline Summer	RFG 0210	5	425
<input type="radio"/> RBOB	RFG 0210	15	500

Summer Spec



The offset values are the negative values of the G-D data column in above D86 report. When we recalculate using these offsets we get:

Quick Reference Guide



As expected, exactly the correct values, because the RFG 0210 offsets are used.

Quick Reference Guide

7.D86 Physical Data Input

Physical D86 data can be entered to determine the setpoints and the offsets of a sample type:

1.0

Still Peaks Run

Offset

Sample	[kPa]
pmp Reformat	29.44

Enter D86 data of pmp Reformat

Input	Volume	T [°C]	StDev
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			

D86 data

Evaporated

Recovered

Cross Check

Save

Type: Reformat
Tsample (5°C) Loss % -
Theater (500°C) Recovery %
Power Residue %
Time before IBP Time rec 5%
Rate (4.5mL/min) Time last 5%

15.09.2004 17:40:18 c:\vhpchem\1\data\b\015F1501.D\DHA_pmp Reformat.csv

The entry of data should be evaporated data beginning with IBP a number of volumes and temperatures ending with FBP. Hit save to store this data with the samplename for later reference. Every following calculation with the same sample name will be compared to these distillatin data. Sample names and data can be removed from the file D86evaporateDHAtoD86.txt. Only evaporated data and crosscheck data can be entered, the recovered data entry is not yet implemented.